



May, 7th 2025

Dear Valued Customer,

RE: BD Vacutainer® PushButton Blood Collection Sets & BD Vacutainer® UltraTouch™ Push Button Blood Collection Sets SKU availability

BD Specimen Management is committed to providing the best possible product experience to our customers. This letter contains valuable information regarding SKU updates on specific BD products to continue to best serve customers.

To enhance efficiency in our manufacturing plants and ensure stable product availability for our customers, we have temporarily paused production of certain blood collection set SKUs. Details are provided in Table 1 below. This change is already in effect and will continue through September 2025. BD will notify you when the paused SKUs become available for order.

Please work with your local BD representative to ensure a smooth transition to the new SKU.

We understand the important role that our products play in ensuring accurate diagnostic test results for your patients as well as safety and confidence for your clinical staff. For more information about these changes, please contact your local BD representative.

Sincerely,

Cormac Smyth
BD SM EMEA Marketing Director

Table 1: SKU Details

	Paused SKU	SKU Description	Alternative SKU*	Alternative SKU Description
1.	367335	25 G x .75 in. BD Vacutainer® Ultra Touch™ Push Button Blood Collection Set with 7 in. tubing and pre-attached holder	367363	25 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and luer adapter
2.	367341	25 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 7 in. tubing and pre-attached holder	367363	25 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and luer adapter
3.	367352	21 G x .75 in. BD Vacutainer® Push Button Blood Collection Set	368689	21 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and pre-attached holder

		with 12 in. tubing and pre-attached holder		
4.	367354	23 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 7 in. tubing and pre-attached holder	368685	23 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 7 in. tubing and pre-attached holder
5.	367355	21 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 7 in. tubing and pre-attached holder	368689	21 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and pre-attached holder
6.	368657	21 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 7 in. tubing and pre-attached holder	368689	21 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and pre-attached holder
7	368658	21 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 12 in. tubing and pre-attached holder	368689	21 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and pre-attached holder
8	368659	25 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 12 in. tubing and pre-attached holder	368687	25 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and pre-attached holder
9	368684	25 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 7 in. tubing and pre-attached holder	368687	25 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and pre-attached holder
10	368686	21 G x .75 in. BD Vacutainer® Ultra Touch™ Push Button Blood Collection Set with 7 in. tubing and pre-attached holder	368689	21 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and pre-attached holder
11	367323	25 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 12 in. tubing without luer	367363	25 G x .75 in. BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with 12 in. tubing and luer adapter
12	367326	21 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 12 in. tubing without luer	367324	23 G x .75 in. BD Vacutainer® Push Button Blood Collection Set with 12 in. tubing without luer

**The potential alternative products are BD's suggestions based on its own comparison of the form, features, and functions of the discontinued/paused products and the potential alternative products. Depending on a customer's circumstances, clinical requirements and intended use of the proposed alternative product, such product may not be a clinically appropriate alternative to the discontinued product, and any such determination shall be made by the customer.*